

# GROWING YOUNG GRACEFULLY WITH STOP™

# stop™

Clinical Skin Renewal Device



**STOP™** the clock and reverse signs of ageing with the revolutionary skin renewal device that offers clinical results by increasing dermal collagen production, resulting in skin tightening and reduced fine lines and wrinkles from the comfort of your own home.


**STOP™** is the perfect choice for those concerned with lined or lax skin wishing to turn back vital signs of ageing without turning to invasive surgery techniques or visiting professional clinics.

**STOP™** is an innovative hand held device powered by **TriPollar™ Radio Frequency (RF) technology**. This 3<sup>rd</sup> generation RF technology has already transformed the professional beauty market with anti-ageing treatments that boost skin renewal producing younger, smoother skin. Now the same proven technology is available in a sleek, palm sized device offering clinical results in the home.


Scientifically developed and clinically proven, **TriPollar™** technology focuses four low power beams deep into the dermis; much like a laser lens focuses its light. Emitting **low power high frequency electrical energy** into the dermis, dermal activity is stimulated to tighten collagen fibres and increase new collagen production. Fibroblasts are stimulated increasing collagen replenishment, resulting in dermal thickening. A thicker dermis will make the skin of a woman in her 40's look and feel like she is in her 30's with visible long lasting results - reduction of fine lines and wrinkles and improved skin texture.

From the initial treatment the results of smoother more radiant skin are noticeable. Over a period of 4-6 weeks, repeating the 15 minute treatment twice per week, long lasting effects are achieved.

Increase of Collagen Fibers in Dermal Layer



before TriPollar™ treatment



after TriPollar™ treatment

Note increase in collagen fibers thickness and quantity after treatment course.

Messon staining pictures of skin biopsy during TriPollar™ clinical study.




Clinical studies showed up to 50% increase in dermal thickening resulting in skin tightening and reduction of fine lines and wrinkles.

\*clinical results by STOP™


## STOP™: technological innovation

Through creating a scientifically-proven product for the home featuring **TriPollar™**, the third generation of Radio Frequency technology, **STOP™** puts the power to control the signs of ageing literally in the hands of the consumer.



## STOP™: the clinical proof

The product was tested by Dr. Ghislaine Beilin (M.D. Paris), specialist in anti-ageing and TriPollar™ technology. The trial, which tested 23 women for the duration of 6 weeks, took place in early 2008. Case studies aged 37-64 used **STOP™** 2-3 times a week and were evaluated by the Fitzpatrick Wrinkle Classification System at the end of their treatments. After their first treatment, objective measurements showed on average an 86% improvement in periorbital wrinkles, and a 76% improvements in perioral wrinkles.



The long term wrinkles reduction was on average 35% for Periorbital wrinkles and 40% for Perioral wrinkles with a maximum wrinkles reduction of 67%.

In vitro clinical test performed by the société de recherche GREDECO in collaboration with researchers from the hôpital Pitié-Salpêtrière, Paris France, on skin biopsies show production of new dermal collagen by more than 40% after **STOP™** treatment.

## STOP™: how to use

The product, which redefines contours on the face and neck area includes **STOP™ Preparation Gel** and **STOP™ After Treatment**.

Simply prepare the face by applying **STOP™ Preparation Gel**; connect the device to a power supply, select from a setting of low, medium or high. Progress with the treatment by ensuring all four electrodes touch the skin. Glide the device in even, circular motions across the face for approximately 15 minutes and a warm, pleasant, massaging sensation will be felt. Following the treatment apply the **STOP™ After Treatment** to enhance the clinical results and soothe the skin.

**Time for the clinical way.  
Time to STOP™.**

### Notes:

**STOP™** in the UK is priced **£250.00** (RRP) and available at **Boots** and **Harvey Nichols** nationwide.

**STOP™** in France is priced **290 EUROS** and available at: Comptoir de l'Homme, Dessange, Darty and selected pharmacies.

**STOP™** is also available on-line at: <https://www.stop-age.com/buy-now.php>

**STOP™** is the latest anti-ageing device from **Ultragen**, leaders in radio frequency technology, and is the first product for the home market offering patent pending **TriPollar™** technology.

**STOP™** should not be used by women who are pregnant, have active skin conditions e.g. Psoriasis or suffer from cardiac disorders or skin cancer.

### About TriPollar™

**TriPollar™ RF** technology was developed by **Pollogen** Ltd in 2006 to offer professional aesthetic systems to medical and aesthetic clinics. The company distributes to over 50 countries globally. **TriPollar™ Radio Frequency** is the integral technology behind **regen™** and **apollo™** machines which provide fat reduction, body contouring, and face firming treatments

**Ultragen** Ltd, a sister company to **Pollogen** Ltd, was established in 2007 to bring **TriPollar™ RF** technology to the consumer with products that offer professional treatments and results in the comfort of the home environment.

For more information on **STOP™** visit [www.stop-age.com](http://www.stop-age.com)

or contact:



**Lim Hock Hin**

Sales Representative

Tel: 63382619

[enquiries@clariancy.com](mailto:enquiries@clariancy.com)

**TriPollar™**  
Technology Inside